

Oasis Eye Gel

Hypromellose USP 0.3% (enriched with Carbomer 980)
Sterile Lubricant Eye Gel

Composition:

Each gram contains Hypromellose USP 3 mg (enriched with Carbomer 980 USNF 0.22%).

Preservative: Sodium Perborate 0.01%

Pharmacology:


Oasis Eye Gel promotes corneal wetting by the stabilization and thickening the precorneal tear film and prolonging the tear film breakdown time.

Indications:

Temporary relief of burning, irritation and discomfort due to dryness of the eye or exposure to wind or sun.

Dosage & administration:

Put 1 or 2 drops in the affected eyes as needed or as advised by the registered physician.

**Contraindications:**


Hypersensitivity to any of the ingredients of this product.

Side effects:

Oasis Eye Gel is very well tolerated by users. However, some cases of burning, stinging, allergic reactions or red eyes after instillation have been reported. Transient blurring of vision after administration has also been reported.

Precautions:

- If eye pain, changes in vision, continued redness or irritation of the eye are noticed, **Oasis Eye Gel** must be discontinued and a doctor should be consulted.
- If **Oasis Eye Gel** colour changes or becomes cloudy, the product should not be used.
- It should never be injected subconjunctivally, nor should be directly introduced into the anterior chamber of the eye.
- It is not intended to be used during surgery.

**Special Guideline**

- ❖ If more than one topical ophthalmic product is to be used, other eye drops should be instilled at least 5 to 10 minutes apart.
- ❖ If using two types of eye medications such as solution and suspension, use the solution first and then use suspension.
- ❖ If using eye drops and an eye ointment, use the eye drops first then apply the ointment at least 10 minutes later.
- ❖ Don't share eye drops with another persons because sharing drops could spread germs and infection.

Use in special group:**Use in pregnancy & lactation:**

There is no experience regarding the safety of **Oasis Eye Gel** in human pregnancy or lactation.

Use in children:

Studies in the pediatric population have not been performed.

Use in elderly patients:

There is no indication that dosage needs to be modified for the elderly.

Storage:

- Store below 30°C, in a dry place.
- Do not use longer than 30 days after first opening.
- Keep out of the reach of children.

Packaging:

Oasis Eye Gel: Each tube contains 10 g Sterile lubricant eye gel and a patient information leaflet.

Manufactured by:

UNIDO
Pharmaceutical Limited
Batpara, Narsingdi, Bangladesh.

PM 2019

